

Population Health Diabetes Education: The Role of Digital Health & Patient-Generated Health Data


Malinda Peoples¹, Janice Macleod¹, Shelley Taylor², Mansur Shomali^{1,3}

¹WellDoc, Inc. Baltimore, MD, ²Meritus Health System, Hagerstown, MD, ³MedStar Union Memorial Hospital, Baltimore, MD
77th Scientific Sessions of the American Diabetes Association, San Diego, CA, USA, June 9 – June 13, 2017

Introduction

Digital health solutions represent a new category of self-management tools that leverage the anytime capability of smartphones to provide contextual real-time coaching to create the opportunity to provide ongoing diabetes self-management support and tracking. This provides robust patient-generated health data (PGHD) to inform the patient-provider dialogue for shared decision-making. However these data are only beneficial when they are fully integrated into the clinical workflow. Here we share findings from a real-world use of a type 2 digital health tool in population health.

Classification of Health Data by Source and Frequency¹



Methods

The Diabetes Population Health Program is part of an Integrated Delivery Network in the Mid-Atlantic with 8 primary care practices, an endocrinology practice and an ADA recognized DSMES program. In this quality improvement project, adults with type 2 diabetes were introduced to a clinically validated², FDA-cleared, digital tool (BlueStar[®], WellDoc[®], Inc., Columbia, MD³) as part of a population health initiative. The lead diabetes educator (population health manager) mentored the care team to provide individualized education, barrier identification, and provider support. Starting in 2014, medical providers have been prescribing the BlueStar product for patients with type 2 diabetes. The software is downloaded onto smart phones with the iOS[™] or Android[™] operating systems. Users interact with the product on the phone, tablet, or on a personal computer. Data is encrypted and uploaded to WellDoc servers. Individual patient reports were sent (in-app) to the providers by the patients based on provider preference or patient-identified need. Population reports were generated at a provider, practice, or system level. For the purpose of this analysis, user data was de-identified in accordance with the Health Insurance Portability and Accountability Act of 1996.

Results

Demographics


Data collection period	July 1, 2016 – June 1, 2017
Number of active users	36
Number of providers	20
Female/Male %	64/36
Age <40 years	16%
Age 41-50	14%
Age 51-60	40%
Age >60	30%
Non-insulin users	48%
Basal insulin users	8%
Basal + bolus insulin users	25%
Medications unknown	19%
Mean baseline A1C	9.7%
Mean baseline BP	132/82
Mean BMI	38.6

Population Data and Opportunities

BlueStar Practice Population Summary Report

July 1, 2016 – June 1, 2017

Persistence of Use



Measure persistence of use to understand population's potential for engagement with a digital tool.

Types of Entries

Entry	% of users with entry type
Blood glucose	69
Medication	69
Exercise	33
Carbs	47
Notes	31

Monitor user engagement and support optimal use; develop engagement strategies for educators, coaches, and team.

Standard of Care Measures

Exam	% of users with entries	% of users with current entries
Dilated eye exam	25	14
Influenza vaccine	19	11
Foot exam	17	11
DSME/S	22	19
Pneumonia vaccine	11	5

Track adherence with SOC metrics. Develop practice-based protocols for achievement.

Blood Glucose Population Data

2,734 total entries	In 25 users
13 hypo events	In 2 users
77 hyper events	In 7 users

Identify users with recurrent extreme BG events; trigger "hypo" or "hyper" and patient outreach protocol.

Hypo event BG <70 mg/dL; hyper event BG >300 mg/dL

Individual User Data and Opportunities

16 patients sent a total of 46 SMART Visit Reports to providers during the observation period

Take Action! BlueStar Clinical and Self-Management Summary

Lab	Problem Value / Date	Med Record Value / Date	Status	Self-Care Behaviors	Plan
A1C	8.2% (07/20/2015)	7.4% (10/14/2015)	Out of Date/Missing	Monitoring: Entered 23.5% of BG in real time*	Encourage entering BG in real time*
LDL	---	---	Out of Date/Missing	Entered a note with a BG % of the time*	Encourage patient to add a note with a BG especially if lower or higher than usual
BP	134/94 mmHg (04/12/2015)	130/80 mmHg (09/03/2015)	Out of Date/Missing	Healthy Eating: Entered carbs 1 times	Encourage patient to keep track of carbs
BMI	---	25.4 by /46.5 (09/03/2015)	Out of Date/Missing	Medication: Highest carb: 50 grams; Lowest carb: 50 grams	Encourage patient to keep track of carbs
Dilated Eye Exam	---	04/10/2015	Overdue	Has recorded meds: Has not set medication reminders	Encourage patient to keep track of meds
Comprehensive Foot Exam	---	05/15/2015	Overdue	Being Active: No data available	Encourage patient to keep active and enter other exercise information in BlueStar
Influenza Vaccine	---	09/24/2015	Overdue	Other: Weak	Encourage patient to keep active and enter other exercise information in BlueStar
Pneumonia Vaccine	---	---	Overdue	Entries: Total Minutes	Encourage patient to keep active and enter other exercise information in BlueStar
Tobacco Usage	Quit Smoker 15 yrs 7 mo	In Target	In Target	No data available	Encourage patient to keep active and enter other exercise information in BlueStar

Healthcare provider reviews BG trends, standards of care, and patient self-management behavior; clinical decision support informs conversation for shared decision-making.

SMART Visit Observations

1. Projected A1C is above target
2. Fasting hyperglycemia suggests increase in basal insulin
3. Post-meal hypoglycemia suggests change in meal-time insulin
4. Self-Management: suggest that patient record carbs more frequently

Conclusions

Digital health tools have recently made available a tremendous amount of PGHD. If summarized and presented properly, these data may facilitate population health management. Diabetes educators care coordinators, and health coaches may be able to deploy targeted, protocol-driven interventions to the patients who need them. As shown above, PGHD may also enhance clinical decision-making for the provider at the time of or in between patient visits. Given the prevalence and complexity of diabetes and the limited time providers have with patients, digital health tools and PGHD will be critical in improving care, supporting practice efficiency, and impacting the cost burden.

References

1. Iyer, A. (2017). *Case Study: The IoT and Big Data in Healthcare Unleashing the Next Generation of Value Creation* in R. Krohn, D. Metcalf, P. Salber (Eds.), *Connected Health: Improving Care, Health, and Efficiency with Wearables and IoT Solutions* (pp. 202-210). Boca Rotan, FL: CRC Press.
2. Quinn C et al. Cluster-Randomized Trial of a Mobile Phone Personalized Behavioral Intervention for Blood Glucose Control Diabetes Care September 2011; vol. 34 no. 9: 1934-1942
3. <https://www.bluestardiabetes.com/> and <http://www.welldoc.com>